


Name: GA-EyeGaze-Bracket

Model Number: #16.1193



Instructions ID: #00.5079_H



Fixation set / Schraubensatz

- (A) 1x EyeGaze Adapter
EyeGaze Adapter
- (B) **not included**
nicht inklusive
- (C) 1x Allen wrench
Innensechskant
- (D) 2x Washer
Scheibe
- (E) 2x M4x6 Bolt
M4x6 Schraube
- (F) 6x M4x10 Bolt
M4x10 Schraube
- (G) 2x M3x6 Bolt
M3x6 Schraube
- (H) 2x M3x10 Bolt
M3x10 Schraube
- (I) 4x Rubber foot
Gummifuß
- (J) 2x Magnet
Magnet
- (K) 2x Tape
Klebeband
- (L) 4x M4x10 Bolt
M4x10 Schraube
- (M) 2x M2,5x6 Bolt
M2,5x6 Schraube

1

-  Mount your "eye-tracker" onto the "EyeGaze-bracket"
-  Montiere deinen "eye-tracker" an den "EyeGaze-Halter"


Device	Fixation set / Schraubensatz	Position	Page/Seite
PCEye Mini	J + K	Magnetic	3
Eye Tracker 4C	J + K	Magnetic	3
MyGaze Power	J + K	Magnetic	3
MyGaze	F	4	4
PCEye Explore	G	2	4
Alea CAM30NT	G	1	4
TM5 Mini	G	3	4
PCEye Go	H	2	4
EyeV Skyle	M	5	4

2

-  Mount the "EyeGaze-bracket" on your "device adapter"
-  Montiere den "EyeGaze-Halter" an deinen "Geräte Adapter"


"SLS"/"Uni-Tablet SLS"

Page/Seite 5




"UDS/UGA"/"Uni-Tablet"

Page/Seite 6



"Notebook"

Page/Seite 7



Tobii Dynavox PCEye Mini

tobii

MyGaze Power

myGaze

Tobii Dynavox EyeTracker 4C






tobii



Fixation set / Schraubensatz



Eye-tracker position






Allea CAM30NT **1** **G**



PCEye Go **2** **H**



PCEye Explore **2** **G**



TM5 Mini **3** **G**








MyGaze **4** **F**







EyeV Skyle **5** **M**



i Fixation set / Schraubensatz **F G H M**










"Spigot Link System (SLS)"





"UDS/UGA System"





"Netbook / Notebook"

